

User Informations / Technical Data Issued: May 2020

NovaProtect™ Face Protection Shield

Product Description

Personal Disposal Face Protection solutions, to reduce infections, transmissions and injuries. Superior, ergonomic fit, lightweight, distortion-free view. Quick use with side protection, high wear comfort, breathable, skin-friendly forehead padding with non-slip, ultrasonic welded elastic rubber band for permanent safety at the workplaces.

Produkteigenschaften

- Optical Class I
- Optimal use with single / reusable safety respirator masks
- Superior, ergonomic fit, can be worn over prescription or protection glasses
- Quick use with side protection, distortion-free view
- High wearing comfort, breathable skin-friendly forehead padding with non-slip rubber band
- Skin-Friendly, Latex-Free, OEKO-TEXT® 100 Standard certified, odorless and suitable for Medical use
- Approved and certified according to Appendix II of the PPE (EU) 2016/425, EN 166:2001 against droplet infections and liquid splashes. No protection against mechanical impacts / parts, shocks, splinters or similar
- Face Shield made from polyester (PET) film, forehead padding polyurethane (PUR) foam, elastic band polyamide (PA), recyclable



Application Areas

Healthcare, laboratories, security, commerce, trade, production, public / private local and long-haul transportations

References / Packaging

- NP.ES.05 : NovaProtect Face Shields of 5 pieces in outer box
- NP.ES.25 : NovaProtect Face Shield of 25 pieces in outer box
- NP.ES.100 : NovaProtect Face Shields of 100 pieces in outer box
- GTIN Nr. : 7629999072961
- Produktbeschreibung in 28 Sprachen

Approvals

The Eye Protection product meets the basic health and safety requirements of PPE Regulations (EU) 2016/425.

The product was tested in the development phase by ECS GmbH – European Certification Service, Hüttfeldstrasse 50, 73430 Aalen, Germany, (Test Center No.: 1883).

The product was tested and CE-approved, according to EN 166:2001 as an Eye Protection product. Does not protect against mechanical influences / parts, impacts, splinters and similar.

Marking / Labeling / Classification

- Face Shield Film : marked / labeled with CE, trademark or without marking / labeling
- Foam : unmarked
- Elastic Rubber Band : unmarked

Novadis GmbH
advanced solutions
Eichtalstrasse 55
CH-8634 Hombrechtikon

Phone
+41 44 810 04 44

Web
www.novadis.ch

Email
info@novadis.ch

Information regarding of the components / materials / weight

Face Shield Film	: from rectified* polyester (PET), *min. extrusion lines, glass-like transparency
Padding Foam	: from polyurethane (PUR) foam, light blue
Elastic Rubber Band	: from polyamide (PA), white
Total Weight	: approx. 30 Gr. / 1.06 oz

Dermatological tolerances

All materials used are skin-friendly, latex-free and odorless. A reaction cannot be ruled out for susceptible / sensitive persons.

Cleaning and Disinfection

The Face Shield is Non-Sterile packaged. Depending on the application, disinfection or cleaning of the Face Shield with 99% ethanol (ethyl alcohol) is recommended.

Storage / Shelf life

In original packaging, without direct sunlight, dry, dust-free and storable at room temperatures from -0° C. – +6 ° C., 32°F – 140°F. The shelf life of the Face Protection Shield is not specified.

Disposal

The Face Protection Shield and can be thermally recycled or disposed of in landfills. Depending on the type and extent of contamination of the Face Shield and their components, the applicable laws and regulations apply in the respective country.

Warnings and Limitations

- Do not reuse the Face Protection Shield after use in the medical environments (risk of contamination). Dispose immediately and safely!
- Protects against drops and splashes of liquids. It reduces the risk of droplet infection. Due to its design, full protection against viruses, bacteria or other biological agents are limited.
- Does not provide protection against mechanical influences / parts, impacts, splinters or similar, as well as against flammable, decomposing or chemical substances (e.g. oils, acids), light sources, gases or similar.
- When wearing Face Protection Shield over normal correction glasses, mechanical influences / parts (e.g. bumps, splinters or similar) can be transmitted and thus pose a risk of injury.
- Do not reuse contaminated, damaged, scratched, heavily worn / soiled Face Protection Shield and dispose immediately. Damages may not be visible to the user. If reused after cleaning or disinfection with ethanol 99% (ethyl alcohol), it is at the discretion and responsibility of the user.

User Conformity

An EU declaration of conformity in accordance with Regulation (EU) 2016/425 on PPE and Directive of -EC 2001/95 on general product safety can be found on the following website and on the packaging (provided with a QR code) for this product: [www.nova-protect.com / downloads](http://www.nova-protect.com/downloads)